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GENERAL MAINTENANCE AND OPERATING INSTRUCTIONS FOR METRIA pH ELECTRODES

During shipment, it is possible for air bubbles to move into the glass bulb. To remove the air, shake down the electrode in the same manner as a clinical thermometer until the glass bulb is filled with solution.

1. PREPARATION FOR USE

All electrodes are shipped with the measuring end covered with a soaker teat. Remove the soaker teat from the electrode, rinse thoroughly with de-ionised water, and gently dry with a non-abrasive tissue. The electrode is ready for use.

Make sure you do not scratch the membrane, as this can cause unstable and inaccurate readings.

2. CALIBRATION

Please consult the pH meter manual to ensure correct operation.

Always make sure the electrode has stabilised in your first calibration solution before starting the procedure. A 2-point calibration is preferred: use pH 7 and pH 4 buffers for acidic samples; pH 7 and pH 10 buffers for alkaline samples. Please contact us for calibration instructions.

Incorrect calibration will cause inaccurate readings.

3. ELECTRODE STORAGE

Store in an upright position and keep the tip wet.

Soaking Solutions:

- pH electrodes should be stored in 4 Molar KCl for both overnight and long-term storage.
- Reference and Double Junction pH electrodes should be stored in the same solution as inside the outer reference chamber.

4. CARE AND CLEANING

Clean after each use with de-ionised water or a solvent that dissolves your sample. Rinse thoroughly with de-ionised water and gently dry with a non-abrasive tissue. Place back in the soaking solution.

Coating of the Membrane:

Slow response and non-reproducible measurements are signs that the electrodes have become coated.

If the glass becomes coated, the time taken to make a measurement will increase (normally, 95% of the final reading should be achieved in less than 10 seconds in pH buffers).

Rinse with iso-propanol or the solvent that dissolves your sample. This should remove the coating and restore response speed.

If the solvent rinse does not restore the response, soak in 0.1 Molar HCl for five minutes. Remove and rinse with water, then place in 0.1 Molar NaOH for five minutes. Remove, rinse again, and soak in pH 4 buffer for 10 minutes before use. Exercise great care when handling chemicals.

Coating and Clogging of the Junction:

Reference electrodes use a liquid junction for electrical contact with the solution being measured. If the junction becomes clogged or coated, the reference becomes erratic. Cleaning with iso-propanol or 0.1 Molar HCl periodically will enhance the electrode's performance. If this does not work, place the tip of the electrode in warm KCl (50 °C) for 5 minutes.

5. FILLING SOLUTIONS

Certain combination electrodes are filled with KCl gel. These do not require filling and have no filling hole on the side of the electrode.

Reference Electrodes:

- R1/Ag: 4 Molar KCl/AgCl
- R2: Upper chamber – 4 Molar KCl/AgCl; Lower chamber – 0.1 Molar KNO₃. The lower chamber can be exchanged for a sample-compatible solution.

To fill the electrode, slide down the plastic ring or remove the fill hole plug to reveal the hole in the side. Using a small syringe, fill the outer part of the electrode (via the hole) with the filling solution detailed above, until the filling solution is just level with the hole. Slide the plastic ring back to cover the hole.

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